



8D HIFU

High Intensity
Focused Ultrasound

HIFU

HIFU, a non-invasive treatment, goes deeper into the skin to stimulate the production of your body's own collagen and elastin, for achieving natural-looking results without surgery or recovery time.


3 = 10 years
Session younger




Theory

What is 8D HIFU?

High Intensity Focused Ultrasound (HIFU) delivers the heat energy to SMAS layers and subcutaneous tissues, which can stimulate and renew the skin's collagen, and therefore improves the skin surface and reduces sagging. This technique can be applied to the face as well as the whole body. And, best of all, it works equally well on all skin tones.




01




Before the treatment, collagen in our skin is weak and disorganized, which makes our skin saggy with fine lines and wrinkles.

02



Micro-focused ultrasonic rays or energies pass through the skin tissue without damaging the cells in the skin layers. It stimulates the collagen in the skin.

03



The skin responds to the stimulations and creates new collagen and elastin simultaneously. Thus, the skin becomes tighter and returns to its youthful state after one month of treatment.

Handpiece&Interface

01 5 Cartridge

The 5 cartridges can penetrate different depth of the face and body precisely. This is a multifunctional and low investment, high return machine.



For Facial Treatment



DS10-1.5
(Optional)



DS7-3.0
(Standard)



DS4-4.5
(Standard)

For Body Slimming



DS3-7.0
(Optional)



DS2-13.0
(Optional)



DS 3-7.0



DS 2-13.0

Standard

Handpiece&Interface

02 Friendly Software



2 handles interface, clear to choose



Professional system, easy to set up parameters



8D system, displays exactly the treatment area

Application

Suitable for what kind of **conditions** ?

Treat saggy skins and aging for areas such as forehead, eyebrows, cheekbones, plump cheeks, cheeks, chin, neck, and body.



Advantage

Quick & short

Only 30 mins per facial treatment

SMAS contraction

Collagen remodeling and elastin fiber contraction

Lasting effect

Good results will last at least 2-3 years

None downtime

Skin only becomes red in the first several hours. It will return to normal afterwards

Instant result

After treatment, you will see 20-30% of immediate results. Best results will be achieved in 3-6 months

Totally non-invasive



Benefits



HIFU can profit up to \$500/session

Treatment site	Average Fee	Procedures	Total Revenue
Face	\$300	2	\$600
Arm	\$200	2	\$400
Abdomen	\$300	2	\$600
Thigh	\$300	2	\$600
Calf	\$200	2	\$400

Testimonials



Dr Justin Boey Medical Director

Beauty Clinic Owner

My clients like the treatment a lot. They enjoy the process. I am very happy to see my clients' skin getting better.



Dr Haus Ranu

Beauty Clinic Owner

Our clinic has been using the 8D HIFU machine for 3 years. It can achieve the same results as Ultherapy. Obvious lifting effects can be observed after only one-time treatment. HIFU helps us attract more and more customers. Love it!



Norman Bethune

Beauty Salon Owner

The results are simply incredible, especially in the area around my cheeks and mouth, where I can see a constant smoothness and disappearance of the fine lines and expression lines. Also, the wrinkles around my eyes and forehead have been substantially reduced. You can see results within 3 weeks with Botox.

Before & After



Technical Specification

Irradiation Type	High Intensity Focused Ultrasound
Frequency	4MHz, 7MHz,
Max Energy	2.0J/m2
LCD Screen	15 Inch True Color Touch LCD ADSS System
Standard Cartridges for single line handle	4MHz: 4.5mm; (cheek/neck)
	7MHz 3mm(forehead/cheek/neck)
Standard Cartridges for 12 lines handle	3MHz 7mm(body) 2MHz 13mm(body)
Opitonal cartridges for single line handle	2MHz 13mm; (body) 3MHz 7mm;(body)
	10MHz 1.5mm(eyes area)
Demensions of Machine	45*30*15cm
N.W./G.W.	20/22KG
Voltage	AC110-220V 50/60Hz